

MX-PRO[®]

bariatric transport



stryker[®]

EMS Equipment

3800 E. Centre Ave.
Portage, MI 49002 U.S.A.
t: 269 329 2100 f: 269 329 2315
toll free: 800 784 4336

www.ems.stryker.com

Standard Features

- Transport at load height capability
- Positive action height adjustment
- Easy-to-use release handle design
- Color-coded controls
- Lap belt extensions
- High visibility powder-coated frame
- Lightweight, aluminum construction
- Scientifically optimized lift bar design
- Lower lifting bar
- Seven height positions
- Integrated bumper system
- Lift-capable safety bar
- Perforated litter surface
- One-hand release breakaway head section with safety bar
- Floor-mounted safety hook
- One-hand release, fold-down side rails
- One-hand release, infinite positioning, pneumatically assisted backrest
- Oversized wheels with sealed caster and wheel bearings
- Reflective labeling
- Sealed flat mattress
- Shock positioning
- Two lap belts and one four-point shoulder restraint

Optional Features

- Heavy duty two or three-stage IV poles (patient right or left)
- Base storage tray
- Nonpocketed head end storage
- Height limit kit¹
- Permanent or removable oxygen bottle holders (foot end)
- Rigid head and foot end push/pull handles
- Single or dual wheel lock
- X-frame guards
- Side lift handles
- Tow package

Specifications

Model Number	6083
Height²	
Position 1	13.5 in (34 cm)
Position 2	24.5 in (62.2 cm)
Position 3	27.5 in (69.9 cm)
Position 4	30 in (76.2 cm)
Position 5	32.5 in (82.6 cm)
Position 6	34.5 in (87.6 cm)
Position 7	37 in (94 cm)
Length	
Standard	80 in (203 cm)
Minimum	62 in (157 cm)
Width	29 in (74 cm)
Weight³	111 lb (50 kg)
Wheels	
Diameter	6 in (15.2 cm)
Width	2 in (5.1 cm)
Articulation	
Backrest	2– 73°
Shock Position	+14°
Maximum Weight Capacity⁴	850/1600 lb (385/725 kg)
Minimum Operator Required	
Occupied Cot	2
Unoccupied Cot	1
Recommended Fastener System	
Floor Mount	Model 6370, 6377 6378
Wall Mount	Model 6371
Recommended Loading Height	Up to 33 in (84 cm)

Warranty

- One year parts, labor and travel or two years parts
- Lifetime on all welds*

Extended warranties available.

*5-year service life.

¹ Height limit kit recommended for load heights less than 28 inches (71 cm). High height kit recommended for load heights more than 32 inches (81 cm).

² Height measured from bottom of mattress at seat section to ground level.

³ Cot is weighed without mattress and restraints.

⁴ 1600 lb capacity at the lowest position.

Stryker reserves the right to change specifications without notice.

The MX-PRO Bariatric Transport is designed to be compatible with competitive cot fastener systems.

Specifications are rounded to the nearest whole number. Conversions are calculated before rounding.

The yellow and black color scheme is a proprietary trademark of the Stryker Corporation.

