


Power-LOAD®

Cot Fastener System



3800 E. Centre Avenue
Portage, MI 49002 USA
t: 269 329 2100
toll free: 800 327 0770
ems.stryker.com

Standard Features

- Lifting arms
- Head end lock LED indicators
- Control panel
- Cot release handles
- Foot end lock LED indicator
- Linear transfer system
- Inductive charging
- Operation guide

Optional Features

- External power cord
- Wheel guide
- Mass casualty floor mount assembly
- Mass casualty wall mount assembly

Specifications

Model Number	6390
Length	
Overall Length	95 in (241 cm)
Minimum Length	89.5 in (228 cm)
Width	24.5 in (62 cm)
Weight	
Total Weight	211.5 lb (96.5 kg)
Floor Plate Assembly	16.5 lb (7.5 kg)
Anchor Assembly	23 lb (10.5 kg)
Transfer Assembly	67 lb (30.5 kg)
Trolley Assembly	105 lb (48 kg)
Maximum Lift Capacity*	700 lb (318 kg)
Minimum Operator Required	
Occupied Cot	2
Unoccupied Cot	1
Recommended Loading Height	22 in to 36 in (56 cm to 91 cm)
Battery	12 VDC, 5 Ah Lead Acid Battery (6390-001-468)

Warranty

- One-year parts, labor, and travel

Extended warranties available.

*7-year service life.

*Maximum weight capacity represents patient weight. Safe working load of 870 lb (395 kg) represents the sum of the cot total weight and patient.

Stryker reserves the right to change specifications without notice.

In-service video included with every order.

The Power-LOAD cot fastener system is designed to conform to the Federal Specification for the Star-of-Life Ambulance KKK-A-1822.

The yellow and black color scheme is a registered trademark of Stryker Corporation.

Meets dynamic crash standards with Power-PRO™ XT cot fastening system (AS/NZS-4535, BS EN-1789 with X-restraints and SAE J3027)

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: **Power-LOAD, Stryker**, . All other trademarks are trademarks of their respective owners or holder.

Planned Certifications

IPX6 IEC-60601-1 BS EN-1789
AS/NZS-4535 KKK-A-1822

Printed in U.S.A.

Copyright © 2017 Stryker

Mkt Lit-539 08 JUN 2011 Rev D.1